



How to measure single wrap strap

INSTRUCTIONS ON HOW TO MEASURE SINGLE WRAP STRAP

Step one: Identify single wrap strap




The holes are located on the inside part of the chair frame as shown in (Pic.1)

Step two: Measure strap width



Place tape measure in the middle of the strap.
If it measures between 1 ¼" to 1 ½", order 1 ½" strap.
If it measures between 1 ¾" to 2", order 2" strap.
See (Pic.2)

Step three: Measure strap length



Place tape measure on the outside part of the frame (left side) and measure across to the outside part of the frame to the right side as shown in (Pic.3)